

## Ontorpaccept Biosimilar - Anti-CD47 mAb - Research Grade

### OVERVIEW

---

**PRODUCT REFERENCE** PX-TA2019

**CLONALITY**

**ISOTYPE** Fusion - [SIRPA (signal regulatory protein alpha, SHPS1, SIRP, SIRPalpha) - IGHG1 (Fc  
(Fragment constant)]

**SOURCE**

### PRODUCT INFORMATION

---

**ORIGIN SPECIES**

**EXPRESSION SYSTEM** XtenCHO

**BUFFER** 0.01M PBS, pH 7.4.

**STORAGE CONDITIONS** Store at 4 °C for short term. Freeze and store at -80 °C for long term.  
Addition of 20 to 40 % glycerol may improve cryoprotection.  
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.  
Aliquot the product according to your needs to avoid freeze/thaw cycles.

**ALIAS/SYNONYMS** anti-CD47, Leukocyte surface antigen CD47, Antigenic surface determinant protein OA3,  
IAP, Integrin-associated protein, Protein MER6, MER6

