



Rilonacept Biosimilar - Anti-Interleukin-1 beta mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2021

CLONALITY

ISOTYPE Fusion - [IL1RAP (interleukin-1 receptor accessory protein) - IL1R1 (type 1 interleukin-1 receptor, interleukin 1 receptor, type I)]₂ - IGHG1 Fc (Fragment constant)

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-Interleukin-1 beta, IL1B, Catabolin, IL1F2, IL-1 beta

