



Tenecteplase Biosimilar - Anti-t-plasminogen activator mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2024

CLONALITY

ISOTYPE 0

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 30 mM histidine pH 5.8, 10% sucrose, 0.02% Tween80

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-t-plasminogen activator, tPA, t-PA, Reteplase, Tissue-type plasminogen activator,
Alteplase, PLAT

