

Altrakinecept Biosimilar - Anti-IL-4 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2025

CLONALITY

ISOTYPE Fusion - [IL4R(Interleukin-4 receptor subunit alpha, IL-4 receptor subunit alpha; IL-4R subunit alpha; IL-4R-alpha; IL-4RA)]

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-IL-4, Binetrakin, B-cell stimulatory factor 1, Pitrakinra, Lymphocyte stimulatory factor 1, BSF-1, IL4, Interleukin-4



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.