



Corifollitropin alfa Biosimilar - Anti-Follitropin receptor mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA2026
CLONALITY	Monoclonal Antibody
ISOTYPE	Fusion - Follicle-stimulating hormone (FSH) - CTP (C-terminal peptide of the β -subunit of human chorionic gonadotropin (hCG))

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4.
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-Follitropin receptor, LGR1, Follicle-stimulating hormone receptor, FSHR, FSH-R

