

## Balugrastim Biosimilar - Anti-G-CSF receptor mAb - Research Grade

### OVERVIEW

---

**PRODUCT REFERENCE** PX-TA2033

**CLONALITY**

**ISOTYPE** Fusion - [ALB (albumin, human serum albumin, HSA) 25-609 - CSF3 (colony stimulating factor 3, GCSF, granulocyte colony-stimulating factor, G-CSF)]

**SOURCE**

### PRODUCT INFORMATION

---

**ORIGIN SPECIES** Homo sapiens

**EXPRESSION SYSTEM** XtenCHO

**BUFFER** 0.01M PBS, pH 7.4.

**STORAGE CONDITIONS** Store at 4 °C for short term. Freeze and store at -80 °C for long term.  
Addition of 20 to 40 % glycerol may improve cryoprotection.  
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.  
Aliquot the product according to your needs to avoid freeze/thaw cycles.

**ALIAS/SYNONYMS** anti-G-CSF receptor, CSF3R, GCSFR, G-CSF-R, CD114, Granulocyte colony-stimulating factor receptor

