



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Blisibimod Biosimilar - Anti-BAFF mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2034

CLONALITY

ISOTYPE Fusion - [peptide 16-mer - peptide 19-mer]₂ - IGHG1 Fc (Fragment constant)

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES Homo sapiens

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-BAFF, ZTNF4, BLYS, Tumor necrosis factor ligand superfamily member 13B,
Dendritic cell-derived TNF-like molecule, TALL-1, TNF- and APOL-related leukocyte
expressed ligand 1, BLYS, B lymphocyte stimulator, TNFSF13B, CD257, TNFSF20,
TALL1, B-cell-activ



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.