



Efavaleukin alfa Biosimilar - Anti-TNFRSF18 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2035

CLONALITY

ISOTYPE Fusion - [IGHG1 Fc (Fragment constant) - [IL2 (interleukin 2, IL-2)]2]

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES Homo sapiens

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-TNFRSF18, Activation-inducible TNFR family receptor, GITR, CD357, Tumor necrosis factor receptor superfamily member 18, AITR, Glucocorticoid-induced TNFR-related protein

