



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Torapsel Biosimilar - Anti-LECAM3 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2037

CLONALITY

ISOTYPE Fusion - [SELPLG (selectin P ligand, selectin P glycoprotein ligand, PSGL-1, CD162) 42-89]2 - IGHG1 Fc (Fragment constant)

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES Homo sapiens

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-LECAM3, Leukocyte-endothelial cell adhesion molecule 3, Platelet activation dependent granule-external membrane protein, Granule membrane protein 140, GMP-140, CD62P, SELP, GRMP, PADGEM, CD62 antigen-like family member P, P-selectin, GMRP



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.