



Taldefgrobep alfa Biosimilar - Anti-GDF-8 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2038

CLONALITY

ISOTYPE Fusion - IGHG1 Fc (Fragment constant) - [peptidyl linker - Alternative to antigen receptors (domain scaffold: FN1 F10, fibronectin domain F10)]₂

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES Homo sapiens

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-GDF-8, GDF8, MSTN, Myostatin, Growth/differentiation factor 8

