



Calpurbatug Biosimilar - Anti-DNA-binding protein HU-alpha mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA2054
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4.
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-DNA-binding protein HU-alpha, HU-2, NS2, hupA, b4000, JW3964, DNA-binding protein HU-beta, HU-1, NS1, hupB, hopD, b0440, JW0430, Integration host factor subunit alpha, IHF-alpha, ihfA, hid, himA, b1712, JW1702, Integration host factor subunit beta, I

