



Mevonlerbart Biosimilar - Anti-Allergen Fdl mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA2069
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-kappa
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Felinized
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4.
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-Allergen Fdl, Major allergen I polypeptide chain 1, AG4, Allergen Cat-1, Allergen Fel d I-A, Allergen Fel dl, Fel d 1-A, CH1, Major allergen I polypeptide chain 2, AG4, Allergen Cat-1, Allergen Fel d I-B, Allergen Fdl, Fel d 1-B, CH2

