



Danvilostomig Biosimilar - Anti-PD1 and CTLA4 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA2076
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa_scFv
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4.
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-PD1, Programmed cell death protein 1, Protein PD-1, hPD-1, PDCD1, CD279, CTLA4, CTLA-4

