



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Perenostobart Biosimilar - Anti-ENTPD1 mAb - Research Grade

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PX-TA2081 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG4-kappa |
| SOURCE | |

PRODUCT INFORMATION

| | |
|---------------------------|---|
| ORIGIN SPECIES | Homo sapiens |
| EXPRESSION SYSTEM | XtenCHO |
| BUFFER | 0.01M PBS, pH 7.4. |
| STORAGE CONDITIONS | Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| ALIAS/SYNONYMS | anti-ENTPD1, Ecto-ATP diphosphohydrolase 1, NTPDase 1, Ecto-ATPase 1, Ecto-apyrase, Ecto-ATPDase 1, Ectonucleoside triphosphate diphosphohydrolase 1, CD39, Lymphoid cell activation antigen |



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.