



Polzastobart Biosimilar - Anti-ILT-4 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA2082
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-kappa
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4.
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-ILT-4, CD85 antigen-like family member D, Monocyte/macrophage immunoglobulin-like receptor 10, MIR-10, ILT4, LIR-2, Immunoglobulin-like transcript 4, MIR10, CD85d, Leukocyte immunoglobulin-like receptor subfamily B member 2, Leukocyte immunoglobulin-

