



Ulviprubart Biosimilar - Anti-CLEC15A mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA2092
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4.
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-CLEC15A, Mast cell function-associated antigen, Killer cell lectin-like receptor subfamily G member 1, C-type lectin domain family 15 member A, MAFA, KLRG1, MAFA-like receptor, ITIM-containing receptor MAFA-L, MAFAL

