



Umizortamig Biosimilar - Anti-TNFRSF9, PDL1, CD3e, EGFRvIII mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA2093
CLONALITY	Monoclonal Antibody
ISOTYPE	immunoglobulin scFv-G1-scFv-scFv_L-kappa
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4.
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-TNFRSF9, 4-1BB, CD137, PDL1, PD-L1, B7H1, CD3e, EGFRvIII

