



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Zigakibart Biosimilar - Anti-TNFSF13 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA2099
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-kappa
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4.
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-TNFSF13, TNF- and APOL-related leukocyte expressed ligand 2, APRIL, CD256, Tumor necrosis factor ligand superfamily member 13, TNF-related death ligand 1, TALL- 2, ZTNF2, A proliferation-inducing ligand, TALL2, TRDL-1



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.