



Luspatercept Biosimilar - Anti-GDF11 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2101

CLONALITY

ISOTYPE Fusion - [ACVR2B (activin A receptor type 2B, ActR-IIB, ActRIIB)]₂ - IGHG1 Fc
(Fragment constant)

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-GDF11, BMP-11, BMP11, Growth/differentiation factor 11, GDF-11, Bone morphogenetic protein 11

