



## Ramatercept Biosimilar - Anti-GDF-8 mAb - Research Grade

### OVERVIEW

---

**PRODUCT REFERENCE** PX-TA2102

**CLONALITY**

**ISOTYPE** Fusion - [ACVR2B (activin A receptor type 2B, ActR-IIB, ActRIIB)]<sub>2</sub> - IGHG1 Fc  
(Fragment constant)

**SOURCE**

### PRODUCT INFORMATION

---

**ORIGIN SPECIES**

**EXPRESSION SYSTEM** XtenCHO

**BUFFER** 0.01M PBS, pH 7.4.

**STORAGE CONDITIONS** Store at 4 °C for short term. Freeze and store at -80 °C for long term.  
Addition of 20 to 40 % glycerol may improve cryoprotection.  
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.  
Aliquot the product according to your needs to avoid freeze/thaw cycles.

**ALIAS/SYNONYMS** anti-GDF-8, GDF8, MSTN, Myostatin, Growth/differentiation factor 8

