



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Sotatercept Biosimilar - Anti-EDF mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2103

CLONALITY

ISOTYPE Fusion - [ACVR2A (activin A receptor type 2A, ActR-II, ActRIIA)]₂ - IGHG1 Fc (Fragment constant)

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-EDF, Erythroid differentiation protein, Activin beta-A chain, INHBA, Inhibin beta A chain



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.