



## Ensovibep Biosimilar - Anti-Spike protein S1 mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA2105
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	DARpin
<b>SOURCE</b>	

### PRODUCT INFORMATION

---

#### ORIGIN SPECIES

**EXPRESSION SYSTEM** XtenCHO

**BUFFER** Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.

**STORAGE CONDITIONS** Store at 4 °C for short term. Freeze and store at -80 °C for long term.  
Addition of 20 to 40 % glycerol may improve cryoprotection.  
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.  
Aliquot the product according to your needs to avoid freeze/thaw cycles.

**ALIAS/SYNONYMS** anti-Spike protein S1, Receptor-Binding Domain, RBD, Spike glycoprotein, S glycoprotein, E2, Peplomer protein, S, D614G variant, B.1.351/Beta, P.1/Gamma, B.1.617.2/Delta, AY.2/Delta Plus, AY.4.2/Delta Plus, Lambda(C.37), Mu (B.1.621), Omicron (B.1.1.529,

